

ClearGuard HD antimicrobial barrier caps for preventing haemodialysis catheter-related bloodstream infections

Information for the public

Published: 13 December 2021

www.nice.org.uk

NICE has said ClearGuard HD antimicrobial barrier caps can be used in the NHS to prevent bloodstream infections in people having haemodialysis with a central venous catheter.

NICE has also said more data should be collected on long-term exposure to chlorhexidine, especially in children.

Central venous catheters have caps on the end. ClearGuard caps contain a rod that is coated in the disinfectant chlorhexidine acetate, which reduces the risk of bacteria from your catheter entering your bloodstream.

ISBN: 978-1-4731-4367-8